



NEONATOLOGY
Vascular Access


premicath

1Fr PICC line
for newborns
under 1kg



Value Life

A 1Fr PICC specially designed for newborns under 1kg



A large choice of introducers to suit all users


24G breakaway needle

The smallest split introducer for small veins



24G I.V. cannula

Basic non-peelable introducer



Catheter with fixation system for the cannula


20G split Microflash

Fast detection of blood return



microsite

The ideal introducer kit for difficult insertion, based on the modified Seldinger technique (only available separately - code: 1147.02)



Premicath without stylet

Codes	Catheter (1Fr in PUR)			Introducer		Other feature	Quantity
	Length cm	Flow rate ml/min (pressure 1 bar)	Prim. vol. ml	Type	Size		
1261.080	8	1.6	0.07	Breakaway needle	24G	-	10/box
1261.152	15	0.9	0.08	I.V. cannula	24G	-	10/box
1261.20	20	0.7	0.09	Breakaway needle	24G	-	10/box
1261.21	20	0.7	0.09	Breakaway needle	24G	Supplied with syringe adaptor	10/box
1261.22	20	0.7	0.09	microflash	20G	-	10/box
1261.201	20	0.7	0.09	Without introducer	-	-	10/box
1261.302	30	0.6	0.11	I.V. cannula	24G	-	10/box


*Premicath with stylet

Codes	Catheter (1Fr in PUR)			Introducer		Other feature	Quantity
	Length cm	Flow rate ml/min (pressure 1 bar)	Prim. vol. ml	Type	Size		
1261.102	10	1.4	0.08	I.V. cannula	24G	-	10/box
1261.153	15	0.9	0.08	I.V. cannula	24G	-	10/box
1261.203	20	0.7	0.09	Breakaway needle	24G	-	10/box
1261.207	20	0.7	0.09	I.V. cannula	24G	-	10/box
1261.206	20	0.7	0.09	Without introducer	-	-	10/box
1261.208	20	0.7	0.09	microflash	20G	-	10/box
1261.210	20	0.7	0.09	Without introducer	-	Supplied with protective polyethylene sheath	10/box
1261.306	30	0.6	0.11	Without introducer	-	-	10/box
1261.307	30	0.6	0.11	I.V. cannula	24G	-	10/box
1261.310	30	0.6	0.11	Without introducer	-	Supplied with protective polyethylene sheath	10/box

Neonatal PICC placement kit

- **Saving time and money:** Reducing preparation time, improving stock management and reducing waste.
- **Improving safety:** Reducing the risk of infection, simplifying procedures and traceability.
- **Instruments dedicated to neonates:** Non adhesive drapes respectful of the fragile skin of the newborns, neonatal tourniquet designed to reduce skin trauma...

1 Preparation and disinfection



Illustrated IFU



2 drape towels
38 x 38 cm



4 measuring tapes



Neonatal tourniquet

2 Positioning of drapes




2 Non adhesive drapes



Fenestrated transparent drape
(40 x 40 cm) with easy peel technology

3 Catheter placement



4 Catheter fixation



2 Tegaderm™
4 x 4cm



Adhesive sutures
6 x 38mm

Code	Product
80199.519	Neonatal PICC Placement kit

VYSET® Vyset is a Vygon activity, which customizes specific sets according to your hospital procedures and your needs.

OBSTETRICS NEONATOLOGY ENTERAL

For further information, please contact: marketingbenelux@vygon.com

The specifications shown in this leaflet are for information only and are not, under any circumstances, of a contractual nature.

SA Vygon NV • Haachtsesteenweg 1650 • Chaussée de Haecht 1650 • 1130 BRUSSELS • BELGIUM
Tel.: +32 (0)2 706 09 50 • www.vygon.be
Vygon Nederland BV • Kerkhofstraat 21 • 5554 HG VALKENSWAARD • THE NETHERLANDS
Tel.: +31 (0)40 208 93 80 • www.vygon.nl
 Vygon Benelux
YY19164

VYCON
www.vygon.be
www.vygon.nl